



Recombinant Arabidopsis thaliana DeSI-like protein At4g17486 (At4g17486)

Product Code	CSB-MP821946DOA
Abbreviation	At4g17486
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93VG8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MWVPTLSSSS CSSDERDESS GEAALTPVYL NVYDLTPVNN YLYWFGIGIF HSGIEAHNLE YCYGAHEYPT SGVYEVEPRN CPGFIFRRSV LLGTTSMSRS DFRSYMEKLS RKYHGDYHL IAKNCNHFTE EVCLQLTGKP IPGWINRLAR VGSFCNCLLP ESIQLTAVSA LPERLEFSDE DESNSEASSV SDEEGSEQHL INVADREIVY LQNKPVRLTR EEIP
Source	Mammalian cell
Target Names	At4g17486
Protein Names	Recommended name: DeSI-like protein At4g17486 EC= 3.4.-.-
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.