



Recombinant Arabidopsis thaliana Protein TIFY 9 (TIFY9)

Product Code	CSB-YP821967DOA
Abbreviation	TIFY9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93ZM9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSKATIELDF LGLEKKQTNN APKPKFQKFL DRRRSFRDIQ GAISKIDPEI IKSLLASTGN NSDSSAKSRS VPSTPREDQP QIPISPVHAS LARSSTELVS GTVPMTIFYN GSVSVFQVSR NKAGEIMKVA NEAASKKDES SMETDLSVIL PTTLRPKLFQ QNLEGDLPIA RRKSLQRFLE KRKERLVSTS PYYPTSA
Source	Yeast
Target Names	TIFY9
Protein Names	Recommended name: Protein TIFY 9 Alternative name(s): Jasmonate ZIM domain-containing protein 10 Protein JASMONATE-ASSOCIATED 1
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.