



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g19525 (At1g19525)

<b>Product Code</b>	CSB-MP821972DOA
<b>Abbreviation</b>	At1g19525
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q940Z1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSQNGIFPDI LTATALVHMY SKSGNFERAT EAFENLKSYG LRPDEKIYEA MILGYVNAGK PKLGERLMKE MQAKELKASE EVYMALLRAY AQMGDANGAA GISSSMQYAS DGPLSFEAYS LFVEAYGKAG QVDKAKSNFD EMRKLGHKPD DKCIANLVRA YKGENSLDKA LRLLLQLEKD GIEIGVITYT VLVDWMANLG LIEEAEQLLV KISQLGEAPP FELQVSLCCM YSGVRNEKKT LQALGVLEAK RDQMGPNEFD KVISALKRGG FEKDARRMYK YMEARKFLPS QRLQMDMVAP PRAFGSGSGR VRRFNS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At1g19525
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At1g19525
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.