



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 10 (Os01g0738400, LOC_Os01g53650)

Product Code	CSB-BP822023OFG
Abbreviation	Os01g0738400, LOC_Os01g53650
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94JI5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAYETSSDHQ LAAAEFLAA LQVHLAGAEA SSPTWGGRCA YDEDFMMYEF KVRRCPRSRA HEWTSCPYAH PGEAARRRDP SHVTYTGEPG PDFRVAARAA CPRGSGCPFA HGTFETWLHP SRYRTRPCRS GMLCARPVCF FAHNDKELRI VGDDAAAATP SPRSPFTTSE DSPPPSPMDM KQIVLAMQQM DARKATRSVA PKTDMLQQEL EEDAPELGWV SDLLM
Source	Baculovirus
Target Names	Os01g0738400
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 10 Short name= OsC3H10
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.