



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At3g06570 (At3g06570)

Product Code	CSB-EP822026DOA
Abbreviation	At3g06570
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94K34
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSSMLSPPVK KRKIVAPQSQ SASASFQSLP DDLILSIVAR VPRLYHRTVS LVCKSFRSLL VSPELYKARS VSGHTESCLY LSIACYPDYR MFTLCRKPDQ TLTTSEEEEEK KKSNGYYLAP VPDPSHPVY FSSLVTVGSD IYNIAGSHAS SNVSILDCRS NTWREAPRLG VELTSVSASV LDRKIFVVG M YADDEEESK NDFFEVLDTE THTWDPQPFN CSETKDKFLN CRTAFIDGKF LVKPWIHRGV VAYNSKESRW EPVQTKMAMS MFNDSYCQIH NVIYLAFDGR IRWYDDALNC WGDVQGLLEL GNIPHGPSCV RLADYRGNIA VFWFRYLPDN DDYKWKMIWC AEIALERTS WEIWGKVLWF DPVLTVPADY EFKALAATV
Source	E.coli
Target Names	At3g06570
Protein Names	Recommended name: F-box/kelch-repeat protein At3g06570
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.