



Recombinant Drosophila melanogaster N-alpha-acetyltransferase 60 (CG18177)

Product Code	CSB-MP822118DLU
Abbreviation	CG18177
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95SX8
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAQFTLYNKH SAPPSESTR VDCEHVPLCS INDVQLRFLV PDDLTEVRQL CQEWFPIDYP LSWYEDITSS TRFFALAAVY NLAIIGLIVA EIKPYRNVNK EVIANMSDSD ELYTRLGFP MQDKGILPDS MGRSADVGYI LSLGVHRSHR RNGIGSLLLD ALMNHLTTAE RHSVKAIFLH TLTTNQPAIF FYEKRRFTLH SFLPYYYNIR GKGKDGFTYV NYINGGHPWP TLLDHIKHYA SMVRHTSSLC AWLAGRVQQV VRWFYHKLLT RFNFIE
Source	Mammalian cell
Target Names	Naa60
Protein Names	Recommended name: N-alpha-acetyltransferase 60 Short name= dNaa60 EC= 2.3.1.48 EC= 2.3.1.88 Alternative name(s): NatF catalytic subunit
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.