




Recombinant Arabidopsis thaliana 2-Cys peroxiredoxin BAS1, chloroplastic (BAS1)

Product Code	CSB-YP822134DOA
Relevance	May be an antioxidant enzyme particularly in the developing shoot and photosynthesizing leaf. Involved in the detoxification of alkyl hydroperoxides with reducing equivalents provided through the thioredoxin system.
Abbreviation	Recombinant Mouse-ear cress BAS1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96291
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	KAQADDLPLVGNKAPDFEAEAVFDQEFIKVKLSDYIGKKYVILFFYPLDFTFVCP TEITAFSDRHSEFEKLNTEVLGVSVDVSVFSLAWVQTDRKSGGLGDLNYPLIS DVTKSISKSFGVLIHDQGIALRGLFIIDKEGVIQHSTINNLGIGRSVDETMRTLQA LQYIQENPDEVCPAGWKPGEKSMKPDPKLSKEYFSAI
Research Area	Others
Source	Yeast
Target Names	BAS1
Protein Names	Thiol-specific antioxidant protein A
Expression Region	66-266aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	24.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.