



Recombinant Avian infectious bronchitis virus Nucleoprotein (N)

Product Code	CSB-YP822153AIAH
Abbreviation	N
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96605
Product Type	Recombinant Protein
Immunogen Species	Avian infectious bronchitis virus (strain V18/91) (IBV)
Purity	≥85% (SDS-PAGE)
Sequence	MSAGKCLKFDS PAPILKLSKN TGSTPPKVG G TGQASWFQSL KEKKRTGTPP TFEGSGVPDN SNVKPQFQHG YWKRQHRYKP GKGG RKP VAD AWYFYTG TG PFGDLKWGDS NDDVVVWKAK GADTSKIGNY GVRDPDKFDQ APLRFTEGGP DNNYRWDFIA LNRGRSRNSS AVTSRENSRP GSRDSSRGRQ RSRVDDDLID RAAKIIMQQQ KNGSRISKQK ANEMAERKYH KRAIAPGKRI DEVFGQRRKG QAPNFGDDKM IEEGVKDGRL TAMLNLVPTP HACLLGSMVT AKLQPDGLHV RFSFETVVKR EDPQFANYSK ICDECVDGVG TRPKDDPTPR SRAASKDRNS APATPKQQRA KKVHKKKEEE SSLTEEEEEV NKQLEYDDDDV TDIPNKIDWG EGAFDDINI
Source	Yeast
Target Names	N
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= NC Short name= Protein N
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.