



Recombinant Human Loss of heterozygosity 12 chromosomal region 1 protein (LOH12CR1)

Product Code	CSB-EP822164HU
Abbreviation	LOH12CR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q969J3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGSEQSSEAE SRPNDLNSSV TPSPAKHRAK MDDIVVVAQG SQASRNVSND PDVIKLQEIP TFQPLLKGLL SGQTSPNAK LEK LDSQQVL QLCLRYQDHL HQCAEAVAFD QNALVKRIKE MDLSVETLFS FMQERQKRYA KYAEQIQKVN EMSAILRRIQ MGIDQTVPLL DRLNSMLPEG ERLEPF SMK P DRELRL
Source	E.coli
Target Names	BORCS5
Protein Names	Recommended name: Loss of heterozygosity 12 chromosomal region 1 protein
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.