



Recombinant Human Leucine-rich repeat-containing protein 39 (LRRC39)

Product Code	CSB-BP822195HU
Abbreviation	LRRC39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96DD0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTENVVCTGA VNAVKEVWEK RIKKLNEDLK REKEFQHKLV RIWEERVSLT KLREKVTRED GRVILKIEKE EWKTLPSSLL KLNQLQEWQL HRTGLLKIPE FIGRFQNLIV LDLSRNTISE IPPGIGLLTR LQELILSYNK IKTVPKELSN CASLEKLELA VNRDICDLPQ ELSNLLKLTH LDLSMNDFTT IPLAVLNMPA LEWLDMGSNK LEQLPDTIER MQNLHTLWLQ RNEITCLPQT ISNMKNLGTL VLSNNKLQDI PVCMEEMANL RFVNFRDNPL KLVSLPPSE GTDEEEEREL FGLQFMHTYI QESRRRADHQ VNGSTTLPIIS INTDG
Source	Baculovirus
Target Names	LRRC39
Protein Names	Recommended name: Leucine-rich repeat-containing protein 39 Alternative name(s): Densin hlg
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.