



Recombinant Human RING finger and SPRY domain-containing protein 1 (RSPRY1)

Product Code	CSB-EP822199HU-B
Abbreviation	RSPRY1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96DX4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>QGLL LTLEEHIAHF LGTGGAATTM GNSCICRDDS GTDDSVDTQQ QQAENSAVPT ADTRSQPRDP VRPPRRGRGP HEP RRKKQNV DGLVLDTLAV IRTLVDNDQE PPYSMITLHE MAETDEGWLD VVQSLIRVIP LEDPLGPAVI TLLLDECPLP TKDALQKLTE ILNLNGEVAC QDSSHPAKHR NTSAVLGCLA EKLGPASIG LLSPGILEYL LQCLKLQSHPTVMLFALIAL EKFAQTSSENK LTISESSISD RLV TLESWAN DPDY LKRQVG FCAQWSLDNL FLKEGRQLTY EKVNLSIRA MLNSNDVSEY LKISPHGLEA RCDASSFESV RCTFCVDAGV WYYEVTVVTS GVMQIGWATR DSKFLNHEGY GIGDDEYS CA YDGC RQLI WY NARSKPHIHP CWKEGDTVGF LLDLNEKQMI FFLNGNQLPP EKQVFSSTVS GFFAAASFMS YQQCEFNF GA KPFKYPPSMK FSTFNDYAFL TAE EK IILPR HRRLALLKQV SIRENCCSLC CDEVADTQLK PCGHS DLCMD CALQLET CPL CRKEIVSRIR QISHIS</p>
Source	E.coli
Target Names	RSPRY1
Protein Names	Recommended name: RING finger and SPRY domain-containing protein 1
Expression Region	17-576
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.