



# Recombinant Human CAP-Gly domain-containing linker protein 3 (CLIP3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP822200HU-B  |
| <b>Abbreviation</b>      | CLIP3   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q96DZ5  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MTKTD PAPMA PPRGEEEEEE EEEDEPVPEA PSPTQERRQK PVVHPSAPAP<br>LPKYAFTFF DPNDPACQEI LFDPQTTIPE LFAIVRQWVP QVQHKIDVIG<br>NEILRRGCHV NDRDGLTDMT LLHYACKAGA HGVGDPAAAV RLSQQLLAG<br>ADVTLRSRWT NMNALHYAAY FDVPDLVRVL LKGARPRVFN STCSDFNHGS<br>ALHIAASSLC LGAACKLEH GANPALNRK GQVPAEVVPD PMDMSLDKAE<br>AALVAKELRT LLEEAVPLSC ALPKVTLPNY DNVPGNLMLS ALGLRLGDRV<br>LLDGQKTGTL RFCGTTEFAS GQWVGVELDE PEGKNDGSVG GVRFYFICPPK<br>QGLFASVSKI SKAVDAPPSS VTSTPRTPRM DFSRVTKGGR REHKGKKKTP<br>SSPSLGSLQQ RDGAKAEVGD QVLVAGQKQG IFRFYGKTDF APGYWYGIEL<br>DQPTGKHDGS VFGVRYFTCP PRHGVFAPAS RIQRIGGSTD SPGDSVGAKK<br>VHQVTMTQPK RTFTTVRTPK DIASENSISR LLFCCWFPWM LRAEMQS |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | CLIP3   |
| <b>Protein Names</b>     | Recommended name: CAP-Gly domain-containing linker protein 3 Alternative name(s): Cytoplasmic linker protein 170-related 59 kDa protein Short name= CLIP-170-related 59 kDa protein Short name= CLIPR-59  |
| <b>Expression Region</b> | 1-547   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.