



Recombinant Human DPH3 homolog (DPH3)

Product Code	CSB-EP822214HU
Abbreviation	DPH3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FX2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAVFHDEVEI EDFQYDEDESE TYFYPCPCGD NFSITKEDLE NGEDVATCPS CSLIKVIYD KDQFVCGETV PAPSANKELV KC
Source	E.coli
Target Names	DPH3
Protein Names	Recommended name: DPH3 homolog Alternative name(s): CSL-type zinc finger-containing protein 2 DelGEF-interacting protein 1 Short name= DelGIP1
Expression Region	1-82
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a CSL zinc finger-containing protein that is required for diphthamide biosynthesis. The encoded protein is necessary for the initial step in the modification of a histidine residue in elongation factor-2 to diphthamide. This modified residue is a target for ADP ribosylation by the bacterial toxins diphtheria toxin and Pseudomonas exotoxin A. Alternative splicing results in multiple transcript variants that encode the same isoform.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.