



Recombinant Human Far upstream element-binding protein 3 (FUBP3)

Product Code	CSB-BP822226HU
Abbreviation	FUBP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96I24
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AELVQQGQSA PVGMKAEGFV DALHRVRQIA AKIDSIPHLN NSTPLVDPSV YGYGVQKRPL DDGVGNQLGA LVHQRTVITE EFKVPDKMVG FIIGRGGEQI SRIQAESGCK IQIASESSGI PERPCVLTGT PESIEQAKRL LGQIVDRCRN GPGFHNIDIDS NSTIQEILIP ASKVGLVIGR GGETIKQLQE RTGVKMVMIQ DGPLPTGADK PLRITGDAFK VQQAREMVLE IIREKDQADF RGVRGDFNSR MGGSIEVSV PRFAVGIVIG RNGEMIKKIQ NDAGVRIQFK PDDGISPERA AQVMGPPDRC QHAAHIISEL ILTAQERDGF GGLAAARGRG RGRGDWSVGA PGGVQEITYT VPADKCGLVI GKGGENIKSI NQQSGAHVEL QRNPPPNNDP NLRRFTIRGV PQQIEVARQL IDEKVGGTNL GAPGAFGQSP FSQPPAPPHQ NTFPPRSSGC FPNMAAKVNG NPHSTPVSGP PAFLTQGWGS TYQAWQQPTQ QVPSQQSQPQ SSQPNYSKAW EDYYKKQSHA ASAAPQASSP PDYTMAWAEY YRQQVAFYGQ TLGQAQAHSQ EQ
Source	Baculovirus
Target Names	FUBP3
Protein Names	Recommended name: Far upstream element-binding protein 3 Short name= FUSE-binding protein 3
Expression Region	2-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.