



# Recombinant Human Far upstream element-binding protein 3 (FUBP3)

<b>Product Code</b>	CSB-MP822226HU
<b>Abbreviation</b>	FUBP3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96I24
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AELVQQGQSA PVGMKAEGFV DALHRVRQIA AKIDSIPHLN NSTPLVDPSV YGYGVQKRPL DDGVGNQLGA LVHQRTVITE EFKVPDKMVG FIIGRGGEQI SRIQAESGCK IQIASESSGI PERPCVLTGT PESIEQAKRL LGQIVDRCRN GPGFHNDIDS NSTIQEILIP ASKVGLVIGR GGETIKQLQE RTGVKMVMIQ DGPLPTGADK PLRITGDAFK VQQAREMVL E IIREKDQADF RGVRGDFNSR MGGGSIEVSV PRFAVGIVIG RRGEMIKKI Q NDAGVRIQFK PDDGISPERA AQVMGPPDRC QHAAHIISEL ILTAQERDGF GGLAAARGRG RGRGDWSVGA PGGVQEITYT VPADKCGGLVI GKGGENIKSI NQQSGAHVEL QRNPPPNSDP NLRRFTIRGV PQQIEVARQL IDEKVGGTNL GAPGAFGQSP FSQPPAPPHQ NTFPPRSSGC FPNMAAKVNG NPHSTPVSGP PAFLTQGWGS TYQAWQQPTQ QVPSQQSQPQ SSQPNYSKAW EDYYKKQSHA ASAAPQASSP PDYTMAWAEY YRQQVAFYGQ TLGQAQAHSQ EQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FUBP3
<b>Protein Names</b>	Recommended name: Far upstream element-binding protein 3 Short name= FUSE-binding protein 3
<b>Expression Region</b>	2-572
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.