



# Recombinant Human F-box/LRR-repeat protein 20 (FBXL20)

<b>Product Code</b>	CSB-MP822228HU
<b>Abbreviation</b>	FBXL20
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96IG2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRRDVNGVTK SRFEMFSNSD EAVINKKLPK ELLLRIFSFL DVVTLCRCAQ VSRAWNVLAL DGSNWQRIDL FDFQRDIEGR VVENISKRCG GFLRKLSLRG CLGVGDNALR TFAQNCRNIE VLNLNGCTKT TDATCTSLSK FCSKLRHLDL ASCTSITNMS LKALSEGCPL LEQLNISWCD QVTKDGIQAL VRGCGGLKAL FLKGCTQLED EALKYIGAHC PELVTLNLQT CLQITDEGLI TICRGCHKLQ SLCASGCSNI TDAILNALGQ NCPRLRILEV ARCSQLTDVG FTTLARNCHE LEKMDLEECV QITDSTLIQL SIHCPRLQVL SLSHCELITD DGIRHLGNGA CAHDQLEVIE LDNCPPLITDA SLEHLKSCHS LERIELYDCQ QITRAGIKRL RTHLPNIKVH AYFAPVTPPP SVGGSRQRFC RCCIIL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FBXL20
<b>Protein Names</b>	Recommended name: F-box/LRR-repeat protein 20 Alternative name(s): F-box and leucine-rich repeat protein 20 F-box/LRR-repeat protein 2-like
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.