



Recombinant Human Terminal nucleotidyltransferase 5A (TENT5A)

Product Code	CSB-MP822229HU
Abbreviation	TENT5A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96IP4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEGEGYFAM SEDELACSPY IPLGGDFGGG DFGGGDFGGG DFGGGGSFGG HCLDYCESPT AHCNVLNWEQ VQRLDGILSE TIPIHGRGNF PTLELQPSLI VKVVRRLAE KRIGVRDVRL NGSAAASHVLH QDSGLGYKDL DLIFCADLRG EGEFQTVKDV VLDCLLDFLP EGVNKEKITP LTLKEAYVQK MVKVCNDSDR WSLISLSNNS GKNVELKFVD SLRRQFEFSV DSFQIKLDSL LLFYECSENP MTETFHPTII GESVYGDFQE AFDHLCNKII ATRNPEEIRG GGLLKYNLL VRGFRPASDE IKTLQRYMCS RFFIDFSDIG EQQRKLESYL QNHFGLEDK KYEYLMTLHG VVNESTVCLM GHERRQTLNL ITMLAIRVLA DQNVIPNVAN VTCYYQPAPY VADANFSNYY IAQVQPVFTC QQQTYSTWLP CN
Source	Mammalian cell
Target Names	TENT5A
Protein Names	Recommended name: Protein FAM46A Alternative name(s): HBV X-transactivated gene 11 protein HBV XAg-transactivated protein 11
Expression Region	1-442aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.