



# Recombinant Human HHIP-like protein 1 (HHIPL1)

<b>Product Code</b>	CSB-BP822236HU
<b>Abbreviation</b>	HHIPL1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96JK4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	H PQCLDFRPPF RPTQPLRLCA QYSDFGCCDE GRDAELTRRF WALASRVDAE EWAACAGYAR DLLCQECSPY AAHLYDAEDP FTPLRTVPGL CQDYCLDMWH KCRGLFRHLS TDQELWALEG NLARFCRYLS LDDTDYCFPY LLVNKNLNSN LGHVVADAKG CLQLCLEEVA NGLRNPVAMV HARDGTHRFF VAEQVGLVWA YLPDRSRLGK PFLNISRVVL TSPWEGDERG FLGIAFHPSF QHNRRLYVYY SVGIRSSEWI RISEFRVSED DENAVDHSSE RIILEVKEPA SNHNGGQLLF GDDGYLYIFT GDGGMAGDPF GTFGNAQNKS ALLGKVLRID VDRKERGLPY GIPPDNPFVG DPAAQPEVYA LGVRNMWRCS FDRGDPSSGT GRGRFCGDV GQNKFEEDV VERGGNYGWR AREGFECYDR SLCANTSLND LLPIFAYPHT VGKSVTGGYV YRGCEYPNLN GLYIFGDFMS GRLMSLQENP GTGQWQYSEI CMGHGQTCEF PGLINNYYPY IISFGEDEAG ELYFMSTGEP SATAPRGVVY KIIDASRRAP PGKCQIQPAQ VKIRSRLIPF VPKEKFIPKT RSTPRPTARA PTRAPRRGRP TAAPPAPTPR PARPTQQPGS RRGGRRRRGR LNSASRAFRD GEVRLVRPAG LSSGSGRVEV FVGGRWGTVC DDSWNISGAA VVCRQLGFAY AVRAVKRAEF GQGGSLPILL DDVRCAGWER NLLCQHNGV GTHNCEHDED AGVVC SHQNP DL
<b>Source</b>	Baculovirus
<b>Target Names</b>	HHIPL1
<b>Protein Names</b>	Recommended name: HHIP-like protein 1
<b>Expression Region</b>	20-782
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.