



Recombinant Human Protein FAM84B (FAM84B)

Product Code	CSB-MP822243HU
Abbreviation	FAM84B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96KN1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGNQVEKLTH LSYKEVPTAD PTGVDRDDGP RIGVSYIFSN DDEDVEPQPP PQGPDGGGLP DGGDGPPPPQ POPYDPRLHE VECSVFYRDE CIYQKSFAPG SAALSTYTPE NLLNKCKPGD LVEFVSQAQY PHWAVYVGNF QVVHLHRLEV INSFLTDASQ GRRGRVVNDL YRYKPLSSSA VVRNALAHVG AKERELSWRN SESFAAWCRY GKREFKIGGE LRIGKQPYRL QIQLSAQRSH TLEFQSLEDL IMEKRRNDQI GRAAVLQELA THLHPAEPEE GDSNVARTTP PPGRPPAPSS EEEDGEAVAH
Source	Mammalian cell
Target Names	FAM84B
Protein Names	Recommended name: Protein FAM84B Alternative name(s): Breast cancer membrane protein 101 Protein NSE2
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.