



# Recombinant Human DNA-directed primase/polymerase protein (PRIMPOL)

<b>Product Code</b>	CSB-MP822255HU
<b>Abbreviation</b>	CCDC111
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96LW4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNRKWEAKLK QIEERASHYE RKPLSSVYRP RLSKPEEPPS IWRLFHRQAQ AFNFVKSCKE DVHVFALECK VGDGQRIYLV TTYAEFWFY KSRKNLLHCY EVIPENAVCK LYFDLEFNKP ANPGADGKKM VALLIEYVCK ALQELYGVNC SAEDVLNLDS STDEKFSQHL IFQLHDVAFK DNIHVGNFRLR KILQPALDLL GSEDDDSAPE TTGHGFPHFS EAPARQGFSF NKMFTKATE ESWTSNSKKL ERLGSAEQSS PDL SFLVKN NMGEKHLFVD LGVYTRNRF RLYKSSKIGK RVALEVTEDN KFFPIQSKDV SDEYQYFLSS LVS NVRFSDT LRILTCEPSQ NKQKGVGYFN SIGTSVETIE GFQCSPYPEV DHFVLSLVNK DGIKGGIRRW NYFFPEELLV YDICKYRWCE NIGRAHKSNN IMILVDLNE VVYQKCHDPV CKAENFKSDC FPLPAEVCLL FLFKEEEEF TDEADETR SN ETQNP HKPSP SRLSTGASAD AVWDNGIDDA YFLEATEDAE LAEAAENSL L SYNSEVDEIP DELIIEVLQE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PRIMPOL
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 111
<b>Expression Region</b>	1-560
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.