



# Recombinant Human Adenosine deaminase domain-containing protein 1 (ADAD1)

<b>Product Code</b>	CSB-BP822259HU
<b>Abbreviation</b>	ADAD1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96M93
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASNNHWFQS SQVPSFAQML KKNLPVQPAT KTITPTGWS SESYGLSKMA SKVTQVTGNF PEPLLSKNLS SISNPVLPPK KIPKEFIMKY KRGEINPVSA LHQFAQMQRV QLDLKETVTT GNVMGPYFAF CAVVDGIQYK TGLGQNKKES RSNAAKLALD ELLQLDEPEP RILETSGPPP FPAEPVVLSE LAYVSKVHYE GRHIQYAKIS QIVKERFNQL ISNRSEYLKY SSSLAAFIIE RAGQHEVVAI GTGEYNYSQD IKPDGRVLHD THAVVTARRS LLRYFYRQLL LFYSKNPAMM EKSIFCTEPT SNLLTLKQNI NICLYMNQLP KGSAQIKSQL RLNPHSISAF EANEELCLHV AVEGKIYLV YCPKDGVNRI SSMSSDKLT RWEVLGVQGA LLSHFIQPVY ISSILIGDGN CSDTRGLEIA IKQRVDDALT SKLPMFYLVN RPHISLVPSA YPLQMNLEYK FLNLSNWAQGD VSLEIVDGLS GKITESSPFK SGMSMASRLC KAAMLSRFNL LAKEAKKELL EAGTYHAAKC MSASYQEAKC KLKSYLQQHG YGSWIVKSPC IEQFNM
<b>Source</b>	Baculovirus
<b>Target Names</b>	ADAD1
<b>Protein Names</b>	Recommended name: Adenosine deaminase domain-containing protein 1 Alternative name(s): Testis nuclear RNA-binding protein
<b>Expression Region</b>	1-576
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.