



Recombinant Human WD repeat-containing protein 92 (WDR92)

Product Code	CSB-EP822265HU-B
Abbreviation	WDR92
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MX6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSAFEKQPQII AHIQKGFNYT VFDCKWVPCS AKFVTMGNFA RGTGVIQLYE IQHGDLKLLR EIEKAKPIKC GTFGATSLQQ RYLATGDFGG NLHIWNLEAP EMPVYSVKGH KEIINAIDGI GGLGIGEGAP EIVTGSRDGT VKVWDPRQKD DPVANMEPVQ GENKRDCWTV AFGNAYNQEE RVVCAGYDNG DIKLFDLRNM ALRWETNIKN GVCSELEFDRK DISMKNLVAT SLEGKFHVFD MRTQHPTKGF ASVSEKAHKS TVWQVRHLPQ NRELFLTAGG AGGLHLWKYE YPIQRSKKDS EGIEMGVAGS VSLLQNVTLS TQPISSLDWS PDKRGLCVCS SFDQTVRVLI VTKLNKI
Source	E.coli
Target Names	WDR92
Protein Names	Recommended name: WD repeat-containing protein 92 Alternative name(s): WD repeat-containing protein Monad
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.