



Recombinant Human FERM domain-containing protein 6 (FRMD6)

Product Code	CSB-MP822270HU
Abbreviation	FRMD6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NE9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNKLNFHNNR VMQDRRSVCI FLPNDESLNI IINVKILCHQ LLVQVCDLLR LKDCHLFGLS VIQNNHEVYM ELSQKLYKYC PKEWKKEASK VRQYEVTWGI DQFGPPMIH FRVQYYVENG RLISDRAARY YYYWHLRKQV LHSQCVLREE AYFLLAAFAL QADLGNFKRN KHYGKYFEPE AYFPSWVVS KRGKDYILKHI PNMHKDQFAL TASEAHLKYI KEAVRLDDVA VHYYRLYKDK REIEASLTG LTMRGIQIFQ NLDEEKQLLY DFPWTNVGKL VFGGKFEIL PDGLPSARKL IYYTGCPMRS RHLLQLLSNS HRLYMNLPV LRHIRKLEEN EEKKQYRESY ISDNLDLMD QLEKRSRAG SSAGSMKHKR LSRHSTASHS SSHTSGIEAD TKPRDTGPED SYSSSAIHRK LKTCSSMTSH GSSHTSGVES GGKDRLEEDL QDDEIEMLV DDPRLDLEQMNE ESLEVSPDMC IYITEDMLMS RKLNGHSGLI VKEIGSSTSS SSETVVKLRG QSTDSLPTI CRKPKTSTDR HSLSLDDIRL YQKDFLRIAG LCQDTAQS YFGCGHELDEE GLYCNSCLAQ QCINIQDAFP VKRTSKYFSL DLTHDEVPEF VV
Source	Mammalian cell
Target Names	FRMD6
Protein Names	Recommended name: FERM domain-containing protein 6 Alternative name(s): Willin
Expression Region	1-622
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.