



Recombinant Human Leukocyte receptor cluster member 8 (LENG8)

Product Code	CSB-MP822284HU
Abbreviation	LENG8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96PV6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AANVGDQRS TDWSSQYSMV AGAGRENGME TPMHENPEWE KARQALASIS KSGAAGGSAK SSSNGPVASA QYVSQAEASA LQQQYYQWY QQYNYAYPYS YYYPMPVPG MDESMSYQAP PQQLPQAQPP QPSNPPHGAH TLNSGPQPGT APATQHSQAG PATGQAYGPH TYTEPAKPKK GQQLWNRMKP APGTGGLKFN IQKRPFVTT QSFGSNAEGQ HSGFGPQPNP EKVNQHSQSS ARGNLSGKPD DWPQDMKEYV ERCFTACESE EDKDRTEKLL KEVLQARLQD GSAYTIDWSR EPLPGLTREP VAESPKKKRW EAASSLHPPR GAGSATRGGG APSQRGTPGA GGAGRARGNS FTKFGNRNVF MKDNSSSSST DSRSRSSRS PTRHFRRSDS HSDSDSSYSG NECHPVGRRN PPPKGRGGRG AHMDRGRGRA QRGKRHD LAP TKRSRKKMAA LECEDPEREL KKQKRAARFQ HGHSRRLRLE PLVLQMSSLE SSGADPDWQE LQIVGTCPDI TKHYLRLTCA PDPSTVRPVA VLKSLCMVK CHWKEKQDYA FACEQMKSIR QDLTVQGIRT EFTVEVYETH ARIALEKGDH EEFNQCQTQL KSLYAENLPG NVGEFTAYRI LYYIFTKNSG DITTELAYLT RELKADPCVA HALALRTAWA LGNYHRFFRL YCHAPCMMSGY LVDFKADRER KVALKAMIKT YVVPSSLLPL LFPSFRLVSS CIKLSPCPLL PLPPVLLIF FSPFPPSLSA FFPWFSTGKT VPLPPSSAM
Source	Mammalian cell
Target Names	LENG8
Protein Names	Recommended name: Leukocyte receptor cluster member 8
Expression Region	2-779
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.