



Recombinant Human Serine/threonine-protein kinase H2 (PSKH2)

Product Code	CSB-EP822291HU
Abbreviation	PSKH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QS6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGCGASRKVV PGPPALAWAK HEGQNGAGVG GAGPGPEAAA QAAQRIQVAR FRAKFDPRVL ARYDIKALIG TGSFSRVVRV EQKTTKKPFA IKVMETRERE GREACVSELS VLRRVSHRYI VQLMEIFETE DQVYVMVELA TGGELFDRLI AQGSFTERDA VRILQMVADG IRYLHALQIT HRNLKPENLL YYHPGEESKI LITDFGLAYS GKKSGDWTMK TLCGTPEYIA PEVLLRKPYP SAVDMWALGV ITYALLSGFL PFDESQTRL YRKILKGKYN YTGEPWPSIS HLAKDFIDKL LILEAGHRMS AGQALDHPWV ITMAAGSSMK NLQRAISRNL MQRASPHSQS PGSAQSSKSH YSHKSRHMWS KRNLRIVESP LSALL
Source	E.coli
Target Names	PSKH2
Protein Names	Recommended name: Serine/threonine-protein kinase H2 EC= 2.7.11.1 Alternative name(s): Protein serine kinase H2 Short name= PSK-H2
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.