



Recombinant Human Proline-rich protein 25 (PRR25)

Product Code	CSB-MP822301HU
Abbreviation	PRR25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S07
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MARTDQKPPC RGGCWGQPGH PNTGGAAHP TYHPMGHRPR TCILLRGDQT TGGQAPSREI SLGPWAAGTH FLAISTTPWG RKTPACISEL PTSSGTAQPL ANAVCEVQTV PGPGLRPQGT PAMRAPSHKG TPPTPNPWGP EQPQNRHKHP KKGVTGGPSP PPPAASRYGQ TPGREPRVQA PGLGPCGRPA SGRLLSLHLE KGDGKGTRQR IPLTDAAVGG DRTDIPSAIA AGPARTPDRH GLPIPGSTPT PMVGSGRLGA PVGRSGGGAS ARSSRPSCAN VLLRADASLG TVLSVLWTGQ LSRGWALLPP GDAGRHLETS VISAGVAAGI WLVEPGEEAQ DPATRRRTAPP RRTASPEPPA PGAPLPACPG RIPGAARFGP RSCPLGSPAV LAVTTGWSHR SV
Source	Mammalian cell
Target Names	PRR25
Protein Names	Recommended name: Proline-rich protein 25
Expression Region	1-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.