



# Recombinant Human Niban-like protein 1 (FAM129B)

<b>Product Code</b>	CSB-EP822308HU-B
<b>Abbreviation</b>	FAM129B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96TA1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGDVLSTHLD DARRQHIAEK TGKILTEFLQ FYEDQYGVAL FNSMRHEIEG TGLPQAQLLW RKVPLDERIV FSGNLFQHQE DSKKWRNRFS LVPHNYGLVL YENKAAYERQ VPPRAVINS A GYKILTSVDQ YLELIGNSLP GTTAKSGSAP ILKCPTQFPL ILWHPYARHY YFCMMTEAEQ DKWQAVLQDC IRHCNNGIPE DSKVEGPAFT DAIRMYRQSK ELYGTWEMLC GNEVQILSNL VMEELGPELK AELGPRLK GK PQERQRQWIQ ISDAVYHMVY EQAKARFEEV LSKVQQVQPA MQAVIRTDMD QIITSKEHLA SKIRAFILPK AEVCVRNHVQ PYIPSILEAL MVPTSQGFTE VRDVFFKEVT DMNLNVINEG GIDKLGEYME KLSRLAYHPL KMQSCYEKME SLRLDGLQQR FVSVSSTSVFK QRAQIHMREQ MDNAVYTFET LLHQELGKGP TKEELCKSIQ RVLERVLK KY DYDSSSVRKR FFREALLQIS IPFLLK LAP TCKSELPRFQ ELIFEDFARF ILVENTYEEV VLQTVMKDIL QAVKEAAVQR KHNLYRDSMV MHNSDPNLHL LAEGAPIDWG EEYSNSGGGG SPSPSTPESA TLSEKRRRAK QVVSVVQDEE VGLPFASPE SPPPASP DG V TEIRGLLAQG LRPESPPPAG PLLNGAPAGE SPQPKAPEA SSPPASPLQH LLPGKAVDLG PPKPSDQETG EQVSSPSSH P ALHTTTEDSA GVQTEF
<b>Source</b>	E.coli
<b>Target Names</b>	FAM129B
<b>Protein Names</b>	Recommended name: Niban-like protein 1 Alternative name(s): Meg-3 Melanoma invasion by ERK Short name= MINERVA Protein FAM129B
<b>Expression Region</b>	1-746
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.