



Recombinant Human Niban-like protein 1 (FAM129B)

Product Code	CSB-EP822308HU
Abbreviation	FAM129B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96TA1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGDVLSTHLD DARRQHIAEK TGKILTEFLQ FYEDQYGVVAL FNSMRHEIEG TGLPQAQLLW RKVPLDERIV FSGNLFQHQE DSKKWRNRFS LVPHNYGLVL YENKAAYERQ VPPRAVINS A GYKILTSVDQ YLELIGNSLP GTTAKSGSAP ILKCPTQFPL ILWHPYARHY YFCMMTEAEQ DKWQAVLQDC IRHCNNGIPE DSKVEGPAFT DAIRMYRQSK ELYGTWEMLC GNEVQILSNL VMEELGPELK AELGPRLK GK PQERQRQWIQ ISDAVYHMVY EQAKARFEEV LSKVQQVQPA MQAVIRTDMD QIITSKEHLA SKIRAFILPK AEVCVRNHVQ PYIPSILEAL MVPTSQGFTE VRDVFFKEVT DMNLNVINEG GIDKLGEYME KLSRLAYHPL KMQSCYEKME SLRLDGLQQR FVSVSSTSVFK QRAQIHMREQ MDNAVYTFET LLHQELGKGP TKEELCKSIQ RVLERVLK KY DYDSSSVRKR FFREALLQIS IPFLLK LK LAP TCKSELPRFQ ELIFEDFARF ILVENTYEEV VLQTVMKDIL QAVKEAAVQR KHNLYRDSMV MHNSDPNLHL LAEGAPIDWG EEYSNSGGGG SPSPSTPESA TLSEKRRRAK QVVS VVQDEE VGLPF EASPE SPPPASP DG V TEIRGLLAQG LRPE SPPPAG PLLNGAPAGE SPQP KAAPEA SSPPASPLQH LLPGKAVDLG PPKPSDQETG EQVSSPSSH P ALHTTTEDSA GVQTEF
Source	E.coli
Target Names	FAM129B
Protein Names	Recommended name: Niban-like protein 1 Alternative name(s): Meg-3 Melanoma invasion by ERK Short name= MINERVA Protein FAM129B
Expression Region	1-746
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.