



Recombinant *Schizosaccharomyces pombe* Uncharacterized mannosyltransferase C32H8.08c (pi016, SPBC32H8.08c)

Product Code	CSB-BP822320SXV
Abbreviation	pi016, SPBC32H8.08c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96WW1
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKRYDTLHVW KLIKLLICKF LLFHDFGTAL FTSLVKASML GFHNWRRTYW LYLKKLRPIN DTYDAPFAIG CKNVAYEASQ YPRMNATFMV LARNSDLGCV LSSMNSIERR FNRHFKYPYV FLNDEPFTTE FKKAVKDATD SSIQFGVLDD ELWNFPKDVD KDMIDESIAE QVGVVYANFP SYHKMCRFFS RNFYKHPLMQ QYEWYWRLEP EVTFCDSISY DPFYYMDKHN KVYGYVIAIK ELAKTVPNLF RYTVAHQKIS NLPTTDLWSF FLDKRYETRI KKLKEEQKDQ GYYVLPEPPL NRIDGQIYNL CHFWSNFEIA RLDFYNSKEY NEYVDALENA GGFWTERWGD APVHSLAVGL LLNRSQVHYF RDLGYQHSTI QHCGQEYGCN CDCPFNIPDY ETKPGSCINE WASVMGGFLD E
Source	Baculovirus
Target Names	omh5
Protein Names	Recommended name: Uncharacterized mannosyltransferase C32H8.08c EC=2.4.1.-
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.