



# Recombinant Arabidopsis thaliana Dof zinc finger protein DOF3.3 (DOF3.3)

<b>Product Code</b>	CSB-MP822590DOA
<b>Abbreviation</b>	DOF3.3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8LFV3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MMMETRDPAI KLFGMKIPFP SVFESAVTVE DDEEDDWSGG DDKSPEKVTP ELSDKNNNNC NDNSFNNSKP ETLDKKEATS TDQIESSDTP EDNQQTTPDG KTLKKPTKIL PCPRCKSMET KFCYNNYNI NQPRHFCKAC QRYWTAGGTM RNVPVGAGRR KNKSSSSHYR HITISEALEA ARLDPGLQAN TRVLSFGLEA QQQHVAAPMT PVMKLQEDQK VSNGARNRFH GLADQRLVAR VENGDDCSSG SSVTTSNNHS VDESRAQSGS VVEAQMNNNN NNNMNGYACI PGVPWPYTWN PAMPPPGFYP PPGYPMPFYP YWTIPMLPPH QSSSPISQKC SNTNSPTLGK HPRDEGSSKK DNETERKQKA GCVLVPKTLR IDDPNEAAKS SIWTTLGIKN EAMCKAGGMF KGFDHKTKMY NNDKAENSPV LSANPAALSR SHNFHEQI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CDF3
<b>Protein Names</b>	Recommended name: Dof zinc finger protein DOF3.3 Short name= AtDOF3.3 Alternative name(s): H-protein promoter-binding factor 2a
<b>Expression Region</b>	1-448
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.