



Recombinant *Oryza sativa* subsp. japonica B3 domain-containing protein Os03g0120900 (Os03g0120900, LOC_Os03g02900)

Product Code	CSB-BP822601OFG
Abbreviation	Os03g0120900, LOC_Os03g02900
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LMR9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEFITPIVRP ASAAAGGGEV QESGGRSLAA VEKEHMFDKV VTPSDVGKLN RLVIPKQHA E KYFPLDAASN EKGLLSFED RTGKPWRFRY SYWNSSQSYV MTKGWSRFVK EKRLDAGDTV SFGRGVGEAA RGRLFIDWRR RPDVVAALQP PTHRFAHHP SSIPFAPWAH HHGHGAAAAA AAAAGARFLL PPSSTPIYDH HRRHAHAVGY DAYAAATSRQ VLFYRPLPPQ QQHHPAVVLE SVPVRMTAGH AEPPSAPSKR VRLFGVNLDC ANSEQDHAGV VGKTAPPPPLP SPPSSSSSSS GKARCSLNL D L
Source	Baculovirus
Target Names	Os03g0120900
Protein Names	Recommended name: B3 domain-containing protein Os03g0120900
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.