



Recombinant Drosophila melanogaster Chitinase-like protein Idgf1 (Idgf1)

Product Code	CSB-EP822660DLU
Abbreviation	Idgf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8MM24
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	ASNLICYYDS NSYLRQGLAK MHTNELDLAL QFCTHLVYGY AGLKSGTLEL FSLNVDLDMF YYKDITALRQ KFPQLKILLS VGGDRDVDEA HPNKYVELLE ANRTAQQNF I DSSMILLKRN GFDGLDLAFQ LPRNKPRKVH GSLGSYWKSF KKLFTGDFVV DPQAEHKSQ FTDLVGNIKN AFRSANLMLS LTVLPNVNST WYFDVPKLHP QFDYINLAAF DFLTPLRNPE EADFTAPIFF QDEQNRLPHL NVEFQINYWL QNHCPGQKLN LGIASYGRAW KLSKGSGLSG APIVHETCGV APGGIIQSA EGLLSWPEIC SKLSQNASAQ YRGELAPLRK VTDLTQKYGN YALRPADDNG DFGVWLSFDD PDFAGIKAVY AKGKGLGGIA LFDLSYDDFR GLCTGQKYPI LRSIKYFMG
Source	E.coli
Target Names	Idgf1
Protein Names	Recommended name: Chitinase-like protein Idgf1 Alternative name(s): Imaginal disk growth factor protein 1
Expression Region	21-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.