



Recombinant Human Tetratricopeptide repeat protein 5 (TTC5)

Product Code	CSB-EP822676HU
Abbreviation	TTC5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N0Z6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMADEEEEEVK PILQKLQELV DQLYSFRDCY FETHSVEDAG RKQQDVQKEM EKTLQQMEEV VGSVQGKAQV LMLTGKALNV TPDYSPKAE E LLSKAVKLEP ELVEAWNQLG EYVWKKGDVA AAHTCFSGAL THCRNKVSLQ NLSMVLRLR TDTEDEHSHH VMDSVRQAKL AVQMDVHDGR SWYILGNSYL SLYFSTGQNP KISQQALSAY AQAQEKVDRKA SSNPDLHLNR ATLHKYEESEY GEALEGFSRA AALDPAWPEP RQREQQLLEF LDRLTSLLES KGKVKTKKLQ SMLGSLRPAH LGPCSDGHYQ SASGQKVTLE LKPLSTLQPG VNSGAVILGK VVFSLTTEEK VPFTFGLVDS DGPCYAVMVY NIVQSWGVI GDSVAIPEPN LRLHRIQHKG KDYSFSSVRV ETPLLLLVNG KPQGSSSQAV ATVASRPQCE
Source	E.coli
Target Names	TTC5
Protein Names	Recommended name: Tetratricopeptide repeat protein 5 Short name= TPR repeat protein 5 Alternative name(s): Stress-responsive activator of p300 Short name= Strap
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.