



Recombinant Human CPX chromosomal region candidate gene 1 protein (CPXCR1)

Product Code	CSB-EP822678HU-B
Abbreviation	CPXCR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N123
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSYPTKEGSD TAGNAHKNSE NEPPNDCSTD IESPSADPNM IYQVETNPIN REPGTATSQE DVVPQAAENS ELETEIQKDQ REEDLKEELL LLQTPIPRKL VSHKPLNDRS RSHSGKVEMK ANNFPINHKT RFRLSTSWRV PFINSHEIRS MILHLLCDRY FSQAAGCQNT MWVKRKYIAC LYHPNSFTHH ERAITFRRPS RVHYRPLTE RMTSGKFCKS TDTKGKCRFR AIVRSVLFVS QIQIESIFNI KGFVDILTYI HTMNVMITNT NNGWKYFCPI CGRLFNTYSE LRQHSCSSSG N
Source	E.coli
Target Names	CPXCR1
Protein Names	Recommended name: CPX chromosomal region candidate gene 1 protein Alternative name(s): Cancer/testis antigen 77 Short name= CT77
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.