



Recombinant Human Leucine-rich repeat-containing protein 71 (LRRC71)

Product Code	CSB-MP822705HU
Abbreviation	LRRC71
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4P6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSSEQSAPGA SPRAPRPGTQ KSSGAVTKKG ERAAKEKPAT VLPPVGEEEP KSPEEYQCSG VLETDFELC TRWGYTDFPK VVNRPRPHPP FVPSASLSEK ATLDDPRLSG SCSLNSLESK YVFFRPTIQV ELEQEDSKSV KEIYIRGWKV EERILGVFSK CLPPLTQLQA INLWKVGLTD KLTTFIELL PLCSSTLRKV SLEGNPLPEQ SYHKLMALDS TIAHLSLRNN NIDDRGAQLL GQALSTLHSC NRTLVSLLNLG FNHIGDEGAG YIADGLRLNR SLLWLSLAHN RIQDKGALKL AEVLRAFELT HTEVVERRRL LLEKGTQERS RSPSSSRHGD SKTDREKSQM VGISNSALVD KTDKTQTMKT PKGLGKKKEK SWELAKKEEK LGSGQSPTQG TPKKEDATKA GKGKVTIPEQ KPSRAKGIKI GSREKRSILL ESELVVEATE VVNPLLEPVE HRDGKVFMPG NKVLLHLNLI RNRITEVGLG GFLATVQYQM QFSKAKSASK GPVGLLWLSL AKNCFAPQCP AYAIQELML PRDPIKAKLR EDEAMAFFP
Source	Mammalian cell
Target Names	LRRC71
Protein Names	Recommended name: Leucine-rich repeat-containing protein 71
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.