



Recombinant Human UPF0565 protein C2orf69 (C2orf69)

Product Code	CSB-EP822731HU-B
Abbreviation	C2orf69
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N8R5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SSCSQA RTMNPGGSGG ARCSLSAEVR RRQCLQLSTV PGADPQRSNE LLLLAAAGEG LERQDLP GDP AKEEPQPPPQ HHVLYFPGDV QNYHEIMTRH PENYQWENWS LENVATILAH RFPNSYIWWI KCSRMHLHKF SCYDNFVKS MFGAPEHNTD FGAFKHL YML LVNAFNLSQN SLSKKSLNVW NKDSIASNCR SSPSHTTNGC QGEKVRTCEK SDESAMSFYP PSLNDASFTL IGFSKGCVV L NQLLFELKEA KKDKNIDAFI KSIRTM YWLD GGHS GGSNTW VTYPEVLKEF AQTGIIVH TH VTPYQVRDPM RSWIGKEHKK FVQILGDLGM QVTSQIHFTK EAPSIENHFR VHEVF
Source	E.coli
Target Names	C2orf69
Protein Names	Recommended name: UPF0565 protein C2orf69
Expression Region	25-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.