



Recombinant Human Acyl-coenzyme A thioesterase 4 (ACOT4)

Product Code	CSB-MP822738HU
Abbreviation	ACOT4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9L9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSATLILEPP GRCCWNEPVR IAVRGLAPEQ RVTLRASLRD EKGALFRAHA RYCADARGEL DLERAPALGG SFAGLEPMGL LWALEPEKPF WRFLKRDVQI PFVVELEVLD GHDPEPGRLL CQAQHERHFL PPGVRRQSVR AGRVRATLFL PPGPGPFPGI IDIFGIGGGL LEYRASLLAG HGFATLALAY YNFEDLPNNM DNISLEYFEE AVCYMLQHPQ VKGPGIGLLG ISLGADICLS MASFLKNVSA TVSINGSGIS GNTAINYKHS SIPPLGYDLR RIKVAFSGLV DIVDIRNALV GGYKNPSMIP IEKAQGPILL IVGQDDHNWR SELYAQTVSE RLQAHGKEKP QIICYPGTGH YIEPPYFPLC PASLHRLLNK HVIWGGEPR HSKAQEDAWK QILAFFCKHL GGTQKTAVPK L
Source	Mammalian cell
Target Names	ACOT4
Protein Names	Recommended name: Acyl-coenzyme A thioesterase 4 Short name= Acyl-CoA thioesterase 4 EC= 3.1.2.2 Alternative name(s): PTE-2b Peroxisomal acyl coenzyme A thioester hydrolase Ib Peroxisomal long-chain acyl-CoA thioesterase Ib
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.