



Recombinant Human Sorting nexin-31 (SNX31)

Product Code	CSB-MP822740HU
Abbreviation	SNX31
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9S9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKMHFCIPVS QQRSDALGGR YVLYSVHLDG FLFCRVRYSQ LHGWNEQLRR VFGNCLPPFP PKYYLAMTTA MADERRDQLE QYLQNVTMPD NVLRSDVFVE FLKLAQLNTF DIATKKAYLD IFLPNEQSIR IEIITSDTAE RVLEVVS HKI GLCRELLGYF GLFLIRFGKE GKLSVVKKLA DFELPYVSLG SSEVENCKVG LRKWYMAPSL DSVLMDCRVA VDLLYMQAIQ DIEKGWAKPT QAQRQKLEAF QKEDSQTKFL ELAREVRHYG YLQLDPCTCD YPESGSGAVL SVGNNEISCC ITLPDSQTQD IVFQMSRVKC WQVTFLGTLT DTDGPQRTL NQNLERFQYS EDSCWQWFVI YTKQAFLLSS CLKKMISEKM VKLAAENTEM QIEVPEQSKS KKYHIQQSQQ KDYSSFLSRK SKIKIAKDDC VFGNIKEEDL
Source	Mammalian cell
Target Names	SNX31
Protein Names	Recommended name: Sorting nexin-31
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.