



# Recombinant Human Probable E3 ubiquitin-protein ligase TRIML1 (TRIML1)

<b>Product Code</b>	CSB-EP822741HU-B
<b>Abbreviation</b>	TRIML1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N9V2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSTADLMENL REELTCFICL DYFSSPVTE CGHSFCLVCL LRSWEEHNTPLSCPECWRTL EGPHFQSNER LGRLASIARQ LRSQVLQSED EQGSYGRMPTAKALSDDDEQ GGSFVAQSH GANRVHLSSE AEEHHREKLQ EILNLLRVRRKEAQAVLTHE KERVKLCQEE TKTCKQVVVS EYMKMHQFLK EEEQLQLQLLEQEEKENMRK LRNNEIKLTQ QIRSLSKMIA QIESSSQSSA FESLEEVRGALERSEPLLLQ CPEATTTLS LCRTGMKEM LRKFSTEITL DPATANAYLVLSEDLKSVKY GGSRQQLPDN PERFDQSATV LGTQIFTSGR HYWEVEVGNKTEWVEGICKD SVSRKGNLPK PPGDLFSLIG LKIGDDYSLW VSSPLKGQHVREPVCVGVF LDYESGHIAF YNGTDESLIY SFPQASFQEA LRPIFSPCLPNEGNTDPLT ICSLNSHV
<b>Source</b>	E.coli
<b>Target Names</b>	TRIML1
<b>Protein Names</b>	Recommended name: Probable E3 ubiquitin-protein ligase TRIML1 EC= 6.3.2.- Alternative name(s): RING finger protein 209 Tripartite motif family-like protein 1
<b>Expression Region</b>	1-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.