



Recombinant Human Probable E3 ubiquitin-protein ligase TRIML1 (TRIML1)

Product Code	CSB-EP822741HU
Abbreviation	TRIML1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9V2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSTADLMENL REELTCFICL DYFSSPVTE CGHSFCLVCL LRSWEEHNTPLSCPECWRTL EGPHFQSNER LGRLASLARQ LRSQVLQSED EQGSYGRMPT TAKALSDDEQ GGSFVAQSH GANRVHLSSE AEEHHREKLQ EILNLLRVRR KEAQAVLTHE KERVKLCQEE TKTCKQVVVS EYMKMHQFLK EEEQLQLQLL EQEEKENMRK LRNNEIKLTQ QIRSLSKMIA QIESSSQSSA FESLEEVRGA LERSEPLLLQ CPEATTTLS LCRTGMKEM LRKFSTEITL DPATANAYLV LSEDLKSVKY GGSRQQLPDN PERFDQSATV LGTQIFTSR HYWEVEVGNK TEWEVGICKD SVSRKGNLPK PPGDLFSLIG LKIGDDYSLW VSSPLKGQHV REPVCKVGVF LDYESGHIAF YNGTDESLIY SFPQASFQEA LRPIFSPCLP NEGNTDPLT ICSLNSHV
Source	E.coli
Target Names	TRIML1
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase TRIML1 EC= 6.3.2.- Alternative name(s): RING finger protein 209 Tripartite motif family-like protein 1
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.