



# Recombinant Human Coiled-coil domain-containing protein 79 (CCDC79)

<b>Product Code</b>	CSB-YP822742HU
<b>Abbreviation</b>	CCDC79
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NA31
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESEDTKKTQ EMKTDLNL LLL ECLKYQMDNA FSQKEALVTI HSICQQNSNA SVYFREIGGL MFVKNLAKSS EHSMVKEAAL YTLGAI AEKN VYCQQTLC TS ELFEDLTWFL SNDSNINLKR MSVYVILVLV SNNRTGQTLV RETGCITVLS RLFRTVISKH ELDLSDKNVF QSYQLWSSVC STLCV CVNNP QNDENQM FCC SLFP HANEWL KNCTTPEIIR PICSFIGLTL ANNTYVQKYF VSVGGLDVLS QVLMQLES DS HETLSSAKLA VVVTKTVDAC IADNPTFGIV LSKYHIVSKL LALLLHESLD SGEKFSIMLT LGHCTEDCEE NQYDLFKNNG LPLMIQALTE SQNEELN KAA TFLVHNCKKI TEKLSLSLGE YPFDENETQQ LKDISVKENN LEEHWRKAKE ILHRIEQLER EGN EEEIQRE NYQDNIS SMN ISIQNTWKHL HADRIGR GSK AEDEDKSHSR QLQSYKSHGV MSKACTNDDQ MKTPLKSANP VHACYRESEQ NKTLYKAKSS CNQNLHEETT FEKNFVSQSS DHVFKHPVHI AKNIKQQLPV TDPFTLCSDI INKEVVSFLA TPSCSEMLTY R CSGCIAVEK SLNSRNFSKL LHSCPYQCDR HKVIVEAEDR YKSELRSLI CNKKILLTPR RRQRLSNEST TPGGIKKRRI RKNFTEEEVN YLFNGVKKMG NHWNSILWSF PFQQGRKA VD LAHKYHKLTK HPTCAAS
<b>Source</b>	Yeast
<b>Target Names</b>	TERB1
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 79
<b>Expression Region</b>	1-727
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.