



Recombinant Human Uncharacterized protein C12orf50 (C12orf50)

Product Code	CSB-EP822743HU-B
Abbreviation	C12orf50
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NA57
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEMQQNCSIS CFWETQPLGC VKISCIFYHS KPRNINGLFL PPSSNITLQK EIQEGIPLQS QSQEPLKPQE NISRPIHHPL VLKTNFEEEE EVDEQNDASS LWTKTPEEIE EKRAIKEMCY KSGEYYRFHT PPDILSSKSM TPTAEKQLEK PLENGSELQE GDSLTVPTKL SQYERQGEIK TSLHGKPKTD IAAFENGGGD CYVPQRVIFL GVDESEALTE EKEITISKCS NTKDNKDSPH PKHSLTTRLV PTTHVLNATE NISMKREDP SSMNDVQPVK KPHFKGVKKR KWIYDEPQNF PNSGMQRAVQ APRPQNKMSY HRNNKNRNAE NASYIHVQRD AVRTVALNAP SRSRPTHGSY NKVHANREPK PNLSPDKYTS TSYNDSAWRK RIPFSKTYSK SEKIYPEPRR NGSK
Source	E.coli
Target Names	C12orf50
Protein Names	Recommended name: Uncharacterized protein C12orf50
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.