



Recombinant Human Kelch repeat and BTB domain-containing protein 3 (KBTBD3)

Product Code	CSB-MP822745HU
Abbreviation	KBTBD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NAB2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDNSYAFNQR STCNIGPSEK KNNFLVSEDH GQKILSVLQN FREQNVFYDF KIIMKDEIIP CHRCVLAACS DFFRAMFEVN MKERDDGSVT ITNLSSKAVK AFLDYAYTGK TKITDDNVEM FFQLSSFLQV SFLSKACSDF LIKSINLVNC LQLLSISDSY GSTSLFDHAL HFVQHHSLL FKSSDFLEMN FGVLQKCLES DELNVPEEEM VLKVLSWTK HNLESRQKYL PHLIEKVRLH QLSEETLQDC LFNEESLLKS TNCFDIIMDA IKCVQGSGGL FPDARPSTTE KYIFIHKTEE NGENQYTCY NIKSDSWKIL PQSHLIDLPG SSLSSYGEKI FLTGGCKGKC CRTVRLHIAE SYHDATDQTW CYCPVKNDFF LVSTMKTPRT MHTSVMALDR LFVIGGKTRG SRDIKSLLDV ESYNPLSKEW ISVSPLPRGI YYPEASTCQN VIYVLGSEVE ITDAFNPSLD CFFKYNATTD QWSELVAEFG QFFHATLIKA VPVNCTLYIC DLSTYKVYSF CPDTCVWKGE GSFECAGFNA GAIGIEDKIY ILGGDYAPDE ITDEVQVYHS NRSEWEEVSP MPRALTEFYC QVIQFNKYRD PWFSNLCA
Source	Mammalian cell
Target Names	KBTBD3
Protein Names	Recommended name: Kelch repeat and BTB domain-containing protein 3 Alternative name(s): BTB and kelch domain-containing protein 3
Expression Region	1-608
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.