



Recombinant Human Draxin (DRAXIN)

Product Code	CSB-YP822750HU
Abbreviation	C1orf187
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NBI3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GALAP GTPARNLPEN HIDLPGPALW TPQASHHRRR GPGKKEWGPG LPSQAQDGAV VTATRQASRL PEAEGLLPEQ SPAGLLQDKD LLLGLALPYP EKENRPPGWE RTRKRSREHK RRRDRLRLHQ GRALVRGPSS LMKKAELSEA QVLDAAMEES STSLAPTMMFF LTTFEAAPAT EESLILPVTS LRPQQAQPRS DGEVMPTLDM ALFDWTDYED LKPDGWPSAK KKEKHRGKLS SDGNETSPAЕ GEPCDHHQDC LPGTCCDLRE HLCTPHNRGL NNKCFDDCMC VEGLRCYAKF HRNRRVTRRK GRCVEPETAN GDQGSFINV
Source	Yeast
Target Names	DRAXIN
Protein Names	Recommended name: Draxin Alternative name(s): Dorsal repulsive axon guidance protein Neucrin
Expression Region	26-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.