



Recombinant Human Protein FAM63B (FAM63B)

Product Code	CSB-EP822752HU-B
Abbreviation	FAM63B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NBR6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESSPELQP LEHGVAAGPA SGTGSSQEGL QETRLAAGDG PGVWAAETSG GNGLGAAAAR RSLPDSASPA GSPEVPGPCS SSAGLDLKDS GLESPAAAEA PLRGQYKVTA SPETAVAGVG HELGTAGDAG ARPDLAGTCQ AELTAAGSEE PSSAGGLSSS CSDPSPPGES PSLDSLESFS NLHSFPSSCE FNSEEGAENR VPEEEEGA AV LPGAVPLCKE EEEGETAQVL AASKERFPGQ SVYHIKWIQW KEENTPIITQ NENGPCPLLA ILNVLLLAWK VKLPPMMEII TAEQLMEYLG DYMLDAKPKE ISEIQLNRYE QNMSDAMAIL HKLQTGLDVN VRFTGVRVFE YTPECIVFDL LDIPLYHGWL VDPQIDDIVK AVGNCSYNQL VEKIISCKQS DNSELVSEGF VAEQFLNNTA TQLTYHGLCE LTSTVQEGEL CVFFRNNHFS TMTKYKGQLY LLVTDQGFLT EEKVVWESLH NVDGDGNFCD SEFHLRPPSD PETVYKGQQD QIDQDYLMAL SLQQEQQSQE INWEQIPEGI SDLELAKKLQ EEEDRRASQY YQEQQEAAAA AAAASTQAQQ GQPAQASPSS GRQSGNSERK RKEPREKDKE KEKEKNSCVI L
Source	E.coli
Target Names	MINDY2
Protein Names	Recommended name: Protein FAM63B
Expression Region	1-621
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.