



# Recombinant Human POC1 centriolar protein homolog A (POC1A)

<b>Product Code</b>	CSB-MP822753HU
<b>Abbreviation</b>	POC1A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NBT0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAAPCAEDPS LERHFKGHRD AVTCVDFSIN TKQLASGSMD SCLMVVHMKP QSRAYRFTGH KDAVTCVNFS PSGHLLASGS RDKTVRIWVP NVKGESTVFR AHTATVRSVH FCSDGQSFVT ASDDKTVKVV ATHRQKFLFS LSQHINWVRC AKFSPDGRLI VSASDDKTVK LWDKSSRECV HSYCEHGGFV TYVDFHPSGT CIAAAGMDNT VKVWDVVRTHR LLQHYQLHSA AVNGLSFHPS GNYLITASSD STLKILDLMG GRLLYTLHGH QGPATTVAFS RTGEYFASGG SDEQVMVWKS NFDIVDHGEV TKVPRPPATL ASSMGNLPEV DFPVPPGRGR SVESVQSQPQ EPVSVPQTLT STLEHIVGQL DVLTQTVSIL EQRLTLTEDK LKQCLENQQL IMQRATP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	POC1A
<b>Protein Names</b>	Recommended name: POC1 centriolar protein homolog A Alternative name(s): WD repeat-containing protein 51A
<b>Expression Region</b>	1-407
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.