



Recombinant Human BTB/POZ domain-containing protein KCTD6 (KCTD6)

Product Code	CSB-BP822756HU
Abbreviation	KCTD6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NC69
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDNGDWGYMM TDPVTLNVGG HLYTTSLTTL TRYPDSSLGA MFGGDFPTAR DPQGNFYFDR DGPLFRYVLN FLRTSELTLP LDFKEFDLLR KEADFYQIEP LIQCLNDPKP LYPMDTFEEV VELSSSTRKLS KYSNPVAVII TQLTITTKVH SLLEGISNYF TKWNKHMMMDT RDCQVSFTFG PCDYHQEVSL RVHLMYITK QGFTIRNTRV HHMSERANEN TVEHNWTFGR LARKTDD
Source	Baculovirus
Target Names	KCTD6
Protein Names	Recommended name: BTB/POZ domain-containing protein KCTD6
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.