



Recombinant Human NFATC2-interacting protein (NFATC2IP)

Product Code	CSB-MP822758HU
Abbreviation	NFATC2IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NCF5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEPVGKRGR WSGGSGAGRG GRGGWGGRRGR RPRAQRSPSR GTLDVVSVDL VTDSDEEILE VATARGAADE VEVEPPEPPG PVASRDNSNS DSEGEDRRPA GPPREPVRRR RRLVLDPGEA PLVPVYSGKV KSSLRLIPDD LSELLKLYPPG DEEEAELADS SGLYHEGSPS PGSPWTKTLR TKDKKEKKKT EFLDLDNSPL SPPSPRTKSR THTRALKKLS EVNKRLQDLR SCLSPKPPQG QEQQGQEDEV VLVEGPTLPE TPRLFPLKIR CRADLVRLPL RMSEPLQSVV DHMATHLGVS PSRILLLFGE TELSPTATPR TLKLGVADII DCVVLTSSPE ATETSQQLQL RVQGKEKHQT LEVLSRDSP LKTLMSHYEE AMGLSGRKLS FFFDGTKLSG RELPADLGME SGDLIEVWG
Source	Mammalian cell
Target Names	NFATC2IP
Protein Names	Recommended name: NFATC2-interacting protein Alternative name(s): 45 kDa NF-AT-interacting protein Short name= 45 kDa NFAT-interacting protein Nuclear factor of activated T-cells, cytoplasmic 2-interacting protein
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.