



Recombinant Human Protein LSM14 homolog A (LSM14A)

Product Code	CSB-MP822763HU
Abbreviation	LSM14A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8ND56
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SGGTPYIGS KISLISKAEI RYEGILYITD TENSTVALAK VRSFGTEDRP TDRPIPPRDE VFEYIIFRGS DIKDLTVCEP PKPQCCLPQD PAIVQSSLGS STSSFQSMGS YGPFGRMPTY SQFSPSSLVG QQFGAVGVAG SSLTSFGTET SNSGTLPQSS AVGSAFTQDT RSLKTQLSQG RSPQLDPLR KSPTMEQAVQ TASAHL PAPA AVGRRSPVST RPLPSASQKA GENQEHRRAE VHKVSRPENE QLRNDNKRQV APGAPSAPRR GRGGHRGGRG RFGIRRDGPM KFEKDFDFES ANAQFNKEEI DREFHNKLLK KEDKLEKQEK PVNGEDKGDS GVDTQNSEGN ADEEDPLGPN CYYDKTKSFF DNISCDDNRE RRPTWAEERR LNAETFGIPL RPNRGRGGYR GRGGLGFRGG RGRGGGRGGT FTAPRGFRGG FRGGRRGGREF ADFEYRKTTA FGP
Source	Mammalian cell
Target Names	LSM14A
Protein Names	Recommended name: Protein LSM14 homolog A Alternative name(s): Protein FAM61A Protein SCD6 homolog Putative alpha-synuclein-binding protein Short name= AlphaSNBP RNA-associated protein 55A Short name= hRAP55 Short nam
Expression Region	2-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.