



Recombinant Human Neuroguidin (NGDN)

Product Code	CSB-EP822772HU-B
Abbreviation	NGDN
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8NEJ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AALGVLESD LPSAVTLLKN LQEQVMVTA QVKSLTQKVQ AGAYPTEKGL SFLEVKDQLL LMYLMDLTHL ILDKASGGSL QGHDAVLRLV EIRTVLEKLR PLDQKLYQI DKLIKTAVTG SLENDPLRF KPHPSNMMSK LSSEDEEEDE AEDDQSEASG KKSVMKGVSKK YVPPRLVPVH YDETEAEREK KRLERAKRRA LSSSVIRELK EQYSDAPEEI RDARHPHVTR QSQEDQHRIN YEESMMVRLS VSKREKGRRK RANVMSSQLH SLTHFSDISA LTGGTVHLDE DQNPIKKRKK IPQKGRKKKG FRRRR
Source	E.coli
Target Names	NGDN
Protein Names	Recommended name: Neuroguidin Alternative name(s): EIF4E-binding protein
Expression Region	2-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>